REGLO-CPF *Digital*

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® ceramic piston pumps are designed for pumping and dispensing applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- The circuit between mains supply and pump has to be earthed.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- The pump must be switched off before taking off the pump-head.
- Do not open or remove the housing while the pump is operating.

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Safety precautions

Media containing particles:

We recommend you to rinse the pump-head after every pumping process to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

Do not pump media containing particles that exceed a diameter of 0.8 mm.

Piston stroke or stroke volume must be set large enough with respect to the particle size!

Piston fractures are not covered by our warranty.

Corrosive liquids:

Corrosive liquids may eventually produce leak paths around the sealing surfaces of the pump-head.

The pump should be inspected for leaks on a regular basis.

Dry running:

Before starting to run the pump we advise you to fill the pump-head with liquid.

Ensure that for priming the pump does not run dry for more than 15 sec.

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Safety precautions

- Only new fuses, according to the specifications stated on page 10 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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Pump drive: 2 years

Pump head:

We advise you to read the enclosed operating manual of the pump-head manufacturer.

Piston fractures are not covered by our warranty (see our instructions for flushing the pump-head on page 3)

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions. If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC SA the warranty terms of the specific manufacturer are valid.

The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

The REGLO-CPF *Digital* is a valveless rotary piston dispensing pump.

Package:

- Pump drive ISM 321
- Power cord
- Operating instructions
- Pump-head as ordered

Head 2–25 µl/stroke

- FMI 009 ceramic piston
- FMI 010 steel piston

Head 5–50 µl/stroke

- FMI 005 ceramic piston

Head 10–100 µl/stroke

- FMI 007 ceramic piston

Please check the package

and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
Rear panel

- 1 RS232 IN (female) for operation controlling via PC
- 2 RS232 OUT (male) for cascading up to 8 pumps

Digital input (TTL level)
 Run/Stop, autostart

Analog output
 speed output: 0–9.0 kHz

- 3 Fuse-holder with voltage selector 115/230V
- 4 Mains socket

⚡ Mains voltage

Mains voltage	Fuse rating
220-240 V _{AC}	2 x 400 mA *)
110-120 V _{AC}	2 x 800 mA *)

*) slow-blow type

⚠ Socket/Power cord

Use exclusively the originally supplied power cord.
 The socket must be earthed
 (protective conductor contact).

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
Voltage selector and
window in the fuse-holder

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Voltage setting 115V / 230V and changing the fuses

- ➔ Switch the pump off,
pull out the mains plug.
- 1 Pull out the fuse-holder by
opening it with a small screw-
driver (size 0).
 - 2 Take out the voltage selector
plate. Turn it and re-insert it
into the fuse-holder so that
the required voltage rating is
facing the window of the
fuse-holder.
 - 3 Insert 2 new fuses
230 V_{AC}: 2 x 400 mA (slow-blow)
115 V_{AC}: 2 x 800 mA (slow-blow)
- ⚠ Use always 2 slow-blow
fuses of the same type
complying with the local
mains voltage.
- 4 Shut the fuse-holder. The
voltage rating is visible in the
window.

⚠ Before starting-up


Check if the voltage setting
visible in the window of the fuse-
holder complies with your local
mains voltage.

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

Operating panel

- 1 On/Off switch
 - 2 Digital LED display
 - 3 Run/Stop
 - 4 Double function:
 - Rotation direction
 - Reset
 - 5 Calibration / Saving a set value
 - 6 Reduce value
 - 7 Increase value
 - 8 MODE key
 - for selecting operating mode
 - Escape key for value setting
- 9 LEDs for active operating mode**
- **PUMP • Flow rate**
Pumping by flow rate (ml/min)
 - **DISP • Volume**
Dispensing by volume (ml)
 - **DISP • Time**
Dispensing by time
 - **PAUSE • Time**
Pause time for dispensing at intervals
 - **STROKE • div(μl)**
volume setting of the piston stroke (450 divisions)

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Control keys

a RUN/STOP

- ➔ Starts and stops the pump
- ➔ Starts, interrupts and continues a dispensing cycle

b MODE

- Changes between operating modes (see p. 12)
- Escape key for value setting

c Rotation direction

Changes the rotation direction. By changing the rotation direction or starting the pump, the rotation direction is indicated in the display with a running hyphen.

Reset key for:

- aborting an interrupted dispensing cycle
- resetting a titration volume to ZERO
- Default key for value resetting, e.g. for resetting a flow rate or volume calibration to its default value (see p. 24)

d CAL

Multi-function key for:

- calibrating the flow rate and dispensing volume
- saving a set value
- max. speed (when pump is running)

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Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Max. stroke volume (µl) of the pump-head
- 4 Reverse piston strokes
➔ only displayed if > 0
- 5 Displays the operating mode last used

Resetting the parameters to the default values

- 1 Power switch »OFF«
- 2 Push the **MODE key** whilst switching the pump on

The following parameters are reset to the following default values:

Mode: PUMP • Flow rate
Rotation direction: clockwise

Flow rate:*
depends on piston stroke volume

Dispensing volume:**
depends on piston stroke volume

Dispensing time:
4.5 sec

Pause time:
2.0 sec

* is based on 720 rpm
** is based on 54 revolutions

Dispensing cycles: 0 (infinite)
Reverse piston strokes: 0
Calibration: Default calibration

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Setting the max. stroke volume

Depending on the mounted pump-head the max. stroke volume is 25, 50 or 100 μl. This volume is factory-set. If, for any reason, this volume must be reset, we recommend you to proceed as follows:

- 1 Switch the pump off
- 2 Maintain pressure on the CAL button and switch the pump on again
- 3 Enter the correct value with the ▲▼ keys (25, 50 or 100 μl depending on the pump-head)
- 4 Press the CAL button again; the entered value is saved.

Zero point calibration

Before using the pump the first time, or whenever necessary, we recommend you to calibrate the 0-point of the flow control ring.

- Start the pump in the mode PUMP • Flow rate with the RUN/STOP key and pump the liquid to be delivered.
- While the pump is in operation, turn the knurled adjustment ring to the right until the pump stops delivering the liquid (the motor is still running!).
- Loosen the fixing screw of the revolving flow control ring and turn it to the value '0'.
- Tighten the fixing screw again.
- Set the required stroke volume (see pages 16 and 17).

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Setting the stroke volume

The flow rate is determined by the stroke volume and the rotation speed.

$$\Rightarrow \text{Flow rate} = \text{Stroke volume} \times \text{speed}$$

This combination offers a great number of possibilities to set the required flow rate.

➔ However, there is no recommendation of how to set the stroke volume and rotation speed in an „optimum“ way.

Basically the following applies:

- **Small stroke volume/ high speed:**
 - low pulsation
 - better resolution of the flow rate and dispensing volume
 - tendency to cavitation (insufficient flow-through of the liquid in the pump-head)
- **Large stroke volume/ low speed:**
 - for elevated viscosity
 - for differential pressures
 - reduced resolution of the flow rate and dispensing volume

For setting the stroke volume see page 17

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

Setting the stroke volume

For setting the stroke volume, you proceed as follows:

- 1 Change to mode STROKE●div(μl)
- 2 Press the CAL button. The stroke volume in μl's starts to blink in the digital LED read-out.
- 3 Enter the required stroke volume in μl's with the ▲▼ key
- 4 Apply again the CAL button. The value to be adjusted with the flow control rings appears in the display read-out (0.0–450.0 div.).

➔ Turn the knurled adjustment ring till the correct value on the fixed and revolving flow control ring equals the digital LED read-out. (The fixing screw must not be loosened!)

- 5 Change to mode PUMP • Flow rate or DISP•Volume

➔ For accurate dispensing we recommend you not to set the stroke volume below 10% of the max. value.

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
Pumping by flow rate

- 1 Change mode to PUMP • Flow rate
- 2 Enter the required flow rate with the ▲▼ keys (displayed in ml/min)
Start pump with RUN/STOP
➔ The flow rate can also be adjusted while the pump is running.
➔ For an accurate flow rate we recommend to calibrate the pump. (see following paragraph)
- 3 For fast filling or emptying the system, maintain pressure on the CAL key (while the pump is running).

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Calibrating the flow rate

- 1 Change the mode to PUMP • Flow rate. Set the required flow rate with the ▲▼ keys
➔ Check: have you entered the required stroke volume in the mode STROKE•div(μl) and adjusted the revolving flow control ring accordingly? (see page 15)
- 2 Change the mode to DISP•Time
– Enter 60 secs. with the ▲▼ keys
– Start the pump with the RUN/STOP-key
➔ The pump stops automatically after 60 secs.
➔ Measure the dispensed liquid by volume or weight

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Calibrating the flow rate

Continuation:

Calibrate the ascertained value as follows:

3 Return to the mode PUMP • Flow rate

4 Press the CAL-button (displayed value blinks)
Enter the measured value (balance) with the ▲▼ keys

5 Save with the CAL-key. The pump changes to the mode STROKE•div(μl) and displays the value to be adjusted with the flow control rings (0.0–450.0 div.)

➔ Turn the knurled adjustment ring till the correct value on the fixed (red) and revolving flow control ring (black) equals the digital LED read-out. A complete revolution equals 100 divisions (max. 450 divisions). The fixing screw must not be loosened!

6 Change the mode to PUMP • Flow rate. The pump is now calibrated and ready.

7 Start the pump with the RUN/STOP-key.

➔ Setting the default calibration: see page 24

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
Titrating


➔ For an accurate titration we recommend you to calibrate the pump (see page 18).

- 1 Switch the pump off.
- 2 Connect the titration handle to RS232 OUT.
- 3 Keep ▲ key pressed and switch power 'ON' (displayed value blinks)
 - ➔ Change to 'titr' with the ▲▼ key
- 4 Save with the CAL-key
- 5 Change the mode to PUMP • Flow rate
 - ➔ '0' appears in the display (titration volume).
- 6 Set the required flow rate with the ▲▼ keys. (is indicated in ml/min)
 - ➔ Afterwards the reading in the display returns to '0'.
- 7 Start the pump with the titration handle and fill the tubing line.
 - ➔ Maintain pressure on the push button.
- 8 Reset the titration volume to '0' by using the ◆ key.
- 9 Start the pump with the titration handle and carry out the titration.
 - ➔ The push button must be pressed down during the titration.
 - ➔ The titrated volume is shown in the display
- 10 The volume can be reset to '0' by using the ◆ key.

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Please note

A dispensing cycle can be interrupted by pressing the RUN/STOP key.

- the remaining dispensing time blinks
- by pressing the RUN/STOP key again, the dispensing cycle is continued
- by pressing the -key an interrupted dispensing cycle can be stopped for good

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Dispensing by time

The dispensing time can be entered from 0.1s to 999h.

- 1 Change the mode to DISP•Time
- 2 Enter the required dispensing time with the keys
➔ The display accelerates when pressure on the or key is maintained.

Seconds: 0.1–899.9 (in 0.1s steps)

Minutes: 15'–899' (in 1min steps)

Hours: 15h–999h (in 1h steps)

- 3 Start pump with RUN/STOP

➔ With the keys the dispensing time can be changed while the dispensing cycle is in process.

➔ During the dispensing process the flow rate (PUMP • Flow rate) can be retrieved by pressing the MODE-key. The flow rate can be changed during the dispensing process by using the keys.

➔ Very short dispensing times can result in non-repeatable dispensing volumes.

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Please note

A dispensing cycle can be interrupted by pressing the RUN/STOP key.

- the volume dispensed up to now is blinking
- by pressing the RUN/STOP key again, the dispensing cycle is continued
- by pressing the key the dispensing cycle is stopped for good

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Dispensing by volume

- 1 Change the mode to DISP • Volume
- 2 Use the keys for entering the required dispensing volume. Two display readings are available:
 Microlitre: e.g. 12.5 μl = **12.5μ**
 Millilitre: e.g. 230.5 ml = **230.5**
 ➔ For dispensing accurately we recommend you to calibrate the pump (see page 23).
- 3 Start pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

- ➔ The dispensing speed can be adjusted in the mode PUMP • Flow rate.
- ➔ With the keys the dispensing volume can also be changed even during the dispensing process. The newly entered set point appears shortly in the display.
- ➔ During the dispensing process the flow rate (PUMP • Flow rate) can be retrieved by pressing the MODE key.
- ➔ The flow rate can be changed even during the dispensing process by using the keys.

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
Calibrating the volume

- Change the MODE to DISP • Volume
➔ Use the ▲▼ keys for entering the required dispensing volume
- Start with RUN/STOP
➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes
➔ Measure the dispensed liquid by volume or weight
- Stay in the same mode (DISP • Volume) and press the CAL-key (display blinks)
- Enter the effectively dispensed volume with the ▲▼ keys
- Save with the CAL button. The pump changes to the mode STROKE • div(μl) and displays the value to be adjusted with the flow control rings (0.0–450.0 div.)
➔ Turn the knurled adjustment ring till the correct value on the fixed (red) and revolving flow control ring (black) equals the digital LED read-out. A complete revolution equals 100 divisions (max. 450 divisions). The fixing screw must not be loosened!
➔ Change to the mode DISP • Volume. The pump is now ready for a calibrated dispensing process.
- Start pump with RUN/STOP
➔ Depending on the application conditions repeated calibration may be necessary.

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Default calibration of flow rate

- 1 Change mode to PUMP • Flow rate.
- 2 Press the CAL button (displayed value blinks)
- 3 Press the key. (Default value blinks).
- 4 Confirm with the CAL button. The pump changes to the mode STROKE•div(μl) and displays the value to be adjusted with the flow control rings (0.0–450.0 div.). The fixing screw must not be loosened!

➔ At the same time, the volume is also reset to the default value.

Default calibration of volume

- 1 Change mode to DISP • Volume
- 2 Press the CAL button (displayed value blinks)
- 3 Press the key. (Default value blinks).
- 4 Confirm with the CAL button The pump changes to the mode STROKE•div(μl) and displays the value to be adjusted with the flow control rings (0.0–450.0 div.). The fixing screw must not be loosened!

➔ At the same time, the flow rate is also reset to the default value.

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
Please note

For entering the number of
dispensing cycles see page 28

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Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time

- 1 Go to mode PAUSE • Time
Enter the pause time
(between 0.1s and 999h)
with ▲▼ keys.
- 2 Press MODE key until
both LEDs for PAUSE and
DISP•Time light up
simultaneously.
➔ Enter the dispensing time
(0.1s to 999h) with ▲▼ keys
- 3 Start with the RUN/STOP key.
➔ The pump dispenses now at
intervals.
➔ Active mode blinks.
- 4 Pushing the RUN/STOP key
again interrupts the
dispensing process or pause
period.
➔ The display blinks.

The user can now proceed as
follows:

- ➔ By pressing the RUN/STOP
key again, the dispensing or
pause period is continued.
- ➔ By pressing the ◀ key the
pump is stopped for good.

- ➔ With the ▲▼ keys the dispen-
sing or pause time can be changed
even during the dispensing process.
The newly entered set point
appears shortly on the display.
➔ During the dispensing process
the flow rate can be retrieved
with the MODE button and
adjusted with the ▲▼ keys.

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
Please note

For entering the number of
dispensing cycles see page 29.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time

- 1 Go to the mode PAUSE • Time.
Enter the pause time
(between 0.1s and 999h)
with ▲▼ keys.
- 2 Press MODE key until
both LEDs for PAUSE and
DISP • Volume light up
simultaneously.
➔ Enter the dispensing
volume with the ▲▼ keys.
- 3 Start with the RUN/STOP key.
➔ The pump dispenses now at
intervals.
➔ Active mode blinks.
- 4 Pushing the RUN/STOP key
again interrupts the dispen-
sing process or pause period.
➔ The display blinks.
The user can now proceed as
follows:
➔ By pressing the RUN/STOP
key again, the dispensing or
pause period is continued.
➔ By pressing the ▶ key the
pump is stopped for good.
➔ With the ▲▼ keys the dispen-
sing volume or pause time can be
changed even during the dispen-
sing process. The newly entered
set point appears shortly on the
display.
➔ During the dispensing process
the flow rate can be retrieved
with the MODE button and
adjusted with the ▲▼ keys.

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PUMP ○ Flow rate
DISP ● Volume
DISP ○ Time
PAUSE ● Time
STROKE ○ div(μl)

Number of dispensing cycles

For intermittent dispensing (by time or volume), the number of dispensing cycles can be entered.

- 1 Power switch »OFF«
- 2 Push the ▼ key whilst switching the pump on
- 3 Enter the number of dispensing cycles with the ▲▼ keys.
➔ 0 = continuous dispensing
- 4 Save with the CAL-key
- 5 Change into one of the two modes for dispensing at intervals
– DISP•Time / PAUSE•Time or
– DISP•Volume / PAUSE•Time
- 6 Start the pump with the RUN/STOP-key

➔ During the pauses the display shows the remaining number of dispensing cycles.

➔ During long pauses, the remaining interval time can be checked with the CAL-key.

➔ At any time, the dispensing process or pause period can be interrupted with the RUN/STOP key.

➔ By pressing the RUN/STOP key again the dispensing process or pause is resumed with the remaining amount of dispensing cycles.

➔ An interrupted dispensing cycle can be stopped for good by using the ♦ key.

➔ The number of dispensing cycles remains saved even when the pump is switched off.

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Drip-free dispensing

Reverse piston strokes allow drip-free dispensing (1 to 100 strokes)

- 1 Switch the pump off
- 2 Maintain pressure on key and switch the pump on again
- 3 Display blinks
Enter the number of reverse piston strokes with the keys
- 4 Save input with the CAL key

Reverse piston strokes are active in the following modes:

- a Dispensing by volume (see page 22)
- b Dispensing a volume within a pre-set time (see page 25)
- c Intermittent dispensing by volume (see page 27)

Reject the first dispensing step. Afterwards, the number of reverse piston strokes are compensated.

Deleting the reverse piston strokes

- Repeat step 1 and 2
- Enter figure 0 with the keys
- Save input by pressing the CAL key

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Pumping against pressure

The REGLO-CPF *Digital* can be used for continuous duty at a differential pressure of max. 6.9 bar.

When the pump is idle

We recommend you to flush the pump-head **after every pumping process** to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

Overcurrent protector

The REGLO-CPF *Digital* drive features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

Whenever this situation occurs, **the pump must be switched off immediately**. Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure or piston seizing).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
Please note

The analog interface is integrated in the port 'RS232 OUT' and can be used for remote control (e.g. by foot-switch or hand dispenser).

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Analogue interface

Pin 1: autostart

If this input is connected to the ground, the pump will automatically start in the operating mode set before the pump was switched off.

Pin 4: +5VDC

For connecting an external control unit (current <100mA).

Pin 5: GND

Point of reference for all signals described.

Pin 6: start

Connection for an external trigger device (hand dispenser, foot-switch, etc.) For operating modes, refer to pin 7.

Pin 7: switch mode

- Input open:
Pump starts or stops when the external trigger device is applied
- Connection to ground:
Pump is running as long as the external trigger device is pressed

Pin 8: busy

Inform about pump status.

- Pumpe in operation
= output +5V
- Pumpe in stand-by
= output on ground (The output resistance is 4.7 kΩ)

Pin 9: speed OUT

A square-wave signal with a frequency proportional to the speed (rpm) is at your disposal: 300 pulses are generated per one rotor revolution.

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
⚠ Please note

When switching on the REGLO-CPF *Digital* pump, the serial interface will only respond after a delay of 3 seconds.

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Serial interface

RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1–8). It is factory-set at 1 and can be changed with @. The address will be stored permanently (even after the pump has been switched off). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed by** a carriage return (ASCII 13). The pump confirms most of the commands with an asterix *. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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
RS232 OUT (male)

Pin 2: RS232 Rx

Is used for connecting additional pumps together with pin 3/5.

Pin 3: RS232 Tx

Is used for connecting additional pumps together with pin 2/5.

Pin 5: GND

Point of reference for all signals described.

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
Each pump requires an individual address.
(see page 32)

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009 .

➔ **Each pump must be allocated its own individual address.**

For this procedure, you first have to connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can you cascade them as shown on the left.

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Please note

In order to reach the max. flow rate, the stroke volume must be set to the max. value (25, 50 or 100 µl) depending on the pump-head model.

If, depending on the application, this stroke volume must be reduced (see page 17–18), we recommend you not to go below 10% of the maximum value.

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Foot switch

Order No. ISM 894

This foot switch serves as a start/stop device.
It is very useful when using the pump as a dispenser for filing tubes, bottles, etc..
Both hands are free for handling the bottles and tubing.

Extension cable 9F / 9M

9 pins, female / 9 pins, male for XC0009 (2 m long)
Order No. AG0013

Adaptor 25F / 9M

25 pins, female / 9 pins, male
Order-No. XC0009

Titration dispenser

in preparation

Maintenance

We recommend the user to flush the pump-head thoroughly after pumping in order to prevent the piston from seizing by remaining particles (especially when media has been pumped that contained salt, proteins or particles).

Otherwise, the REGLO-CPF *Digital* requires no particular maintenance or mechanical adjustments.

Avoid pumping media containing abrasive solids.

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
⚠ Please note

For maintenance, service, repairs and warranty of the pump-head we also refer to the enclosed instruction manual of the pump-head manufacturer. In case of doubt that manual shall prevail.

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Service and repairs

Like all pumps, the REGLO-CPF *Digital* contains components that will wear over a period of time.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packing or a packet of equal quality.

Replacement parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask from your ISMATEC® distributor:

■ Piston-/cylinder group	
Type	Order No
H-00CKC-LF	FMI 311
H-00SKY-LF	FMI 310
H-0CKC	FMI 307
H-1CKC	FMI 308

■ Lip seals (2 pieces are needed)	
Type	Order No
H408-00J (RH00)	FMI 325
H408-0A (RH0)	FMI 312
H408-1A (RH1)	FMI 319

- Parts lists
- Replacement parts for pump drive
- Wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

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
Ask for the latest
ISMATEC catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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A Unit of IDEX Corporation
Laboratoriumstechnik GmbH